## Stop using the flotation technique and start weighing salbutamol pressurised metered-dose inhalers without dose counters

Ermindo R. Di Paolo<sup>1</sup>, Stéphanie Spaggiari<sup>1,2</sup>, André Pannatier<sup>1,2</sup>, Farshid Sadeghipour<sup>1,2</sup>, Isabelle Rochat<sup>3</sup>, Gaudenz M. Hafen<sup>3</sup>

- 1 Department of Pharmacy, Lausanne University Hospital, Lausanne, Switzerland
- 2 School of pharmaceutical sciences, University of Geneva, University of Lausanne, Geneva, Switzerland
- 3 Respiratory Unit, Department of Paediatrics, Lausanne University Hospital, Lausanne, Switzerland

## Abstract

**Objectives:** Salbutamol pressurised metered-dose inhalers (pMDIs) are not equipped with dose counters outside the USA. The aim of this study was to describe a simple reproducible method for determining the number of doses remaining in a pMDI based on scale weight.

**Methods:** Sixty salbutamol pMDIs (GlaxoSmithKline) containing 200 doses were bought in six European countries (10 pMDIs, with two source batches per country), including France (Ventoline®), Germany (Sultanol®), Italy, Spain, Switzerland and the United Kingdom (Ventolin®). A laboratory scale (Mettler Toledo XP204®) was used to weigh the canisters. Weights were also obtained with a digital kitchen scale used in our hospital (Ohaus CS 5000®) and one marketed for home use (IKEA Ordning®).

**Results:** With a laboratory scale, the mean weight of the canisters was  $28.61 \pm 0.10$  g after priming and  $14.84 \pm 0.23$  g after 200 puffs. Similar results were obtained with two common digital scales. There was a linear correlation between the mean weights of the canisters measured on the three scales (r = 0.998, p <0.0001). Based on the results, a practical table was produced to determine whether pMDIs are empty or close-to empty.

**Conclusions:** As long as manufacturers neglect to equip all pMDIs with a dose counter, we recommend that the number of doses remaining in salbutamol pMDIs being checked by weighing the pMDIs on a common digital scale. Our results indicate that the canister should be replaced when its weight falls to  $\leq$ 15 g.

Published in: Swiss Med Wkly 2015 doi: 10.4414/smw.2015.14162

Contact: ermindo.di-paolo@chuv.ch